



## MobiMed Release Notes 4.2.43



Software Part Number: [490212288](#)

Software Build Number: Release-4.2.43.0\_19.12.05.2

### Changes specific for SWASFT

- EprExpander missing trigger (Rationale field in ECS) (10218)

Under Cardiac Arrest > Resuscitation > Resuscitation tab:



### Changes specific for NIAS

- Corpus 3 Integration - Transfer data from Corpus (9800)
- Corpus 3 Integration - Import vital values (9801)


- Identify and get access to HCN API (9803)
- Implement HCN lookup simulator (10076)
- Support sharing ePR in a personal device set up - Signing (10130)
- Support sharing ePR in a personal device set up - Dependent variables handling (Time values/Expanders etc) (10146)
- Make installation path configurable (10176)

## General improvements


- ePR
  - Add possibility for FASS button in medication control (10164)

The screenshot shows a medical software interface. On the left, there's a sidebar with various medication names like 'Inj. Adrenalin vid anafylaxi', 'Inj. Adrenalin vid obstruktivitet', and 'Inj. Adrenalin vid hjärtstopp'. Below this is a search bar with a placeholder 'Sök på läkemedel, substans, ATC eller vnr'. To the right of the search bar are buttons for 'Allmänhet', 'Vårdpersonal', and 'Djurläkemedel'. Further down, there's a detailed product page for 'Aerius' (MSD) showing its packaging (two brown tablets) and some text. On the far right, there are sections for 'Receptbelagd' (with a note about symbols) and 'Kontakt' (with a 'Fråga om läk' button). A red arrow points from the 'Tid' button in the search bar to the 'FASS' button on the left.


- Add time button for "At scene" in date control (10159)



- PDet\_patientLookup not cleared on CWS when the user changes patient ID (10213)
- Allowing longer integer/numeric fields with leading +/- (10187)
- VSM
  - Reorder and rename review tabs (10161)



- Adapt monitoring views to FZ-M1's resolution (10162)
- Make it possible to set a default value in admintool to force an ECG cable set (10203)
- Make it possible to view the start and end time in the waveform for an alarm (10204)



- Investigate how to run MobiMed PU as a standard program (10157)
- Update server installation instruction (10170)
- Stop all services before update with Smart Server Installer (10184)
- Make it possible to use the WiFi profiles that Windows handles (10191)
- SOS IntegrationSettings document: Add description of DispatcherServiceName setting (10207)

## Fixes

- VSM:
  - The ECG running waveform stops for a long time if "Add test data" is selected (10219)
  - Monitorview disappears after changing patientype and switching between epr and monitoring (10153)
  - Correct deadlock in arrhythmia control/event navigator (10210)
  - SCR: Selecting "Emergency" reason should show visible label (10234)
  - OSDK with SV as designersetting does not affect personnummer (10126)
  - AdminTool: Hospital Integration view does not support sending same

record to multiple backends (10152)

- SOS Dispatch: Dialog for selecting mission should be updated automatically if new mission detected (10183)
- Correct vcredist validation bug found in Gloucestershire (10205)
- HospitalIntegration startup, should read encrypted settings correctly (10209)
- AdminTool: Start screen, re-add missing info on installed version (10215)
- Correct FZ-M1 screen size when docked/undocked (10225)
- Paratus - CAD service will need two restarts to work when installed on a PU (10237)

*\*Numbers in brackets are ID-numbers from Ortivus' work item tracking tool. Numbers starting with "INC" refers to items in ServiceNow*

## Supported Platforms

<b>Hardware Unit</b>	<b>Operating System</b>
Panasonic FZ-M1	Windows 10 (1709)
Panasonic FZ-G1	Windows 10 (1607, 1709, 1809)
Panasonic H2	Windows 7, Windows 8.0
Panasonic CF-19	Windows 7
Getac RX10H Extreme	Windows 10 (1709, 1809)
Getac RX10H	Windows 10 (1511)

All versions of measuring unit M531 are supported with MobiMed 4.2. Measuring unit 300C is not supported.

The Clinical Workstation (CWS) platform should meet or beat the minimum hardware specification and use a supported operating system.

## Release Information

This is a commercial release on the following markets:

- Norway
- Sweden
- United Kingdom

This release replaces all previous 4.2 releases.

Upgrade to this release can be done from all previous 4.2 releases.

## Release Schedule

Release	Release date	Last feature request
4.2.41	30 Sep 2019	19 Aug 2019
4.2.42	04 Nov 2019	23 Sep 2019
4.2.43	09 Dec 2019	28 Oct 2019
4.2.44	03 Feb 2020	02 Dec 2019
4.2.45	16 Mar 2020	27 Jan 2020
4.2.46	20 Apr 2020	09 Mar 2020
4.2.47	25 May 2020	13 Apr 2020
4.2.48	29 Jun 2020	18 May 2020
4.2.49	31 Aug 2020	22 Jun 2020
4.2.50	05 Oct 2020	24 Aug 2020
4.2.51	09 Nov 2020	28 Sep 2020



### Ortivus

Svärdsvägen 19, Danderyd  
Sweden

[info@ortivus.com](mailto:info@ortivus.com)

You received this email because you signed up on our website or made a purchase from us.

[Unsubscribe](#)

mailerlite