

MobiMed Release Notes 4.49



Software Article Number: [490212424](#)

Software Build Number: Release-4.49_20.08.27.1

Release Updates

The following new features are implemented in 4.49:

1. Updating one patient ePR from multiple MobiMed units at the same time.
2. FHIR-API
3. Change Requests implemented
4. Minor changes
5. Fixes



Contact your Customer Representative for more information.

Updating one patient ePR from multiple MobiMed units at the same time

When working with a patient, it is now possible for two paramedics to share, view, and edit one patient's records using two MobiMed units. This allows e.g. one paramedic to log vital parameters at the same time as a second paramedic log general condition or medicines.

FHIR-API

MobiMed is now able to publish ePR-records using the FHIR-API. Records from a Care Episode including Patient data, Triage, Observations, and the Final Handover report is possible to publish.



MobiMed is supporting v4 of FHIR.

CCN SC302 Respiratory Template

The template developed for SCAS will guide crews through the relevant essential fields, with the stated aims being ease of use for crews and an improvement in cycle time for those incidents. This CCN covers the second template, for patients with respiratory conditions.

CCN SC303 Chest Pain Template

These templates will guide crews through the relevant essential fields, with the stated aims being ease of use for crews and an improvement in cycle time for those incidents. This CCN covers the Chest Pain Template.

Minor Changes and General Improvements

CCN SC304 Covid 19 to presenting condition[11264]

Presenting condition and provisional diagnosis of 'Positive COVID-19'. Please allocate to the 'Respiratory' category.

SWASFT TSW2.0 Cardiac Arrest Page enhancements[11230]

SWASFT's cardiac arrest pages update with new drug shortcut buttons and registration of heart rhythms.

Smart installer updates[10918]

SmartInstaller now defaults to save waveforms to disk. It is also possible to select the location (directory) for waveforms.

Forms for acute missions in Prehospital care for Örebro in Sweden [11020]

Created a patient record for Örebro primary care's mobile teams and reports for them.

Minor change to Remote Update [11172]

Monitoring services will only be installed with monitoring option enabled.

SSIS XML export is now compatible with SQL 2019 in Data Warehouse[11269]

SSIS package is now recompiled to support SQL 2019.

Fixes

- Resolved an issue with the smart server was old DLL's was not cleared[10995]
- Resolved an issue were Incorrect calculation of the present number of unsigned ePRs were not displayed correctly during "hand-over"[11048]
- Resolved an issue were Smart server installer did never set the Ortivus\Update\UpdatePath registry key[11164]
- Resolved an issue in SCAS CCN SC305 Safeguarding Deletion Also Deleting Main Patient Details Screen ID:0101149 [11254]
- Resolved an issue in SCAS Falls Template: Allergies & Medications Deletion. Selecting No Known Allergies hid the fields Important medication" input field and the "Patient's medication Taken to Hospital".Issue ID:0101146 [11255]
- Resolved an issue were An error message was written in the event log during startup if the WaveFileDBPath was set to an empty string in registry. This ECG waves will be stored in the database hence error did affect the function. [11294]

Supported Platforms

Hardware Unit

Getac RX10H

Getac RX10H Extreme

Getac UX10

Operating System

Windows 10 (1511)

Windows 10 (1709, 1809)

Windows 10 (1809)

Lenovo 20L3000	Windows 10 (1809)
Lenovo Thinkpad L490	Windows 10 (1909)
Panasonic CF-19	Windows 7
Panasonic FZ-G1	Windows 10 (1607, 1709, 1809)
Panasonic H2	Windows 7, Windows 8.0
Panasonic FZ-M1	Windows 10 (1709)

All versions of measuring unit M531 are supported with MobiMed 4.2. Measuring unit 300C is not supported.

The Clinical Workstation (CWS) platform should meet or beat the minimum hardware specification and use a supported operating system.

Installation and deployment

For installing this release. Please contact our Service and Support at the below addresses.

- **UK market:** support@ortivus.co.uk
- **Other markets:** support@ortivus.com

Release Information

This is a commercial release on the following markets:

- Norway
- Sweden
- United Kingdom

Upgrade to this release can be done from release 4.2 and onwards.

Please contact your customer representative for more information.

Release Schedule

MobiMed Release Plan

Release	Release date	Last feature request
4.2.44	03 Feb 2020	02 Dec 2019
4.2.45	16 Mar 2020	27 Jan 2020
4.2.46	20 Apr 2020	09 Mar 2020
4.47	25 May 2020	13 Apr 2020
4.48	29 Jun 2020	18 May 2020
4.49	31 Aug 2020	22 Jun 2020
4.50	05 Oct 2020	24 Aug 2020
4.51	09 Nov 2020	28 Sep 2020
4.52	21 Dec 2020	02 Nov 2020
4.53	15 Feb 2021	14 Dec 2020

Ortivus

Svärdsvägen 19, Danderyd
Sweden

+46 8 446 4500
info@ortivus.com

You received this email because you signed up on our
website or made a purchase from us.

[Unsubscribe](#)