

A white square icon containing a black EKG (heart rate) line.

# MobiMed

**Saving vital seconds**



**Compact, robust and intuitive defibrillator**



**CPR assistance**



**360 Joule**



**Proven technology**



**MobiMed connect, communication throughout the care chain**

With over 30 years' experience in cardiology and in close collaboration with our users, Ortivus has developed MobiMed, the leading solution for patient monitoring and record keeping in prehospital care. With MobiMed, ambulance personnel can monitor and treat the patient out in the field, while patient data and vital parameters are shared in real time throughout the care chain. MobiMed life - a part of MobiMed, is a stand-alone defibrillation solution that can be tailored to the organisation's needs.

**We work together with our users to continuously develop MobiMed, for the healthcare of today and tomorrow.**



## MobiMed Life - Heart safe solutions

Ortivus delivers modular, user-friendly and ergonomic solutions focused on quality and patient safety. MobiMed Life is designed for advanced monitoring and resuscitation. It is a modular solution that can be adapted to the organisation's needs. Together with MobiMed Monitor and MobiMed ePR, it constitutes a unique total solution for prehospital healthcare.

### **MobiMed ePR**

MobiMed ePR provides the ambulance paramedic with a structured workflow and support that improves the clinical decision-making process.

### **MobiMed Monitor**

With a total weight of 2.5 kg, the MobiMed monitor can easily be carried out into the field. MobiMed monitors and measures the patient's vital parameters such as ECG, blood pressure and SpO2. Data from MobiMed is automatically transferred to the medical record, which reduces the administrative burden enabling the team to focus on the management of the patient.

### **MobiMed Life**

A defibrillator is not always needed, but when it is, it is crucial that it is close at hand. Data from the defibrillator is automatically transferred to MobiMed monitor and ePR.



## CPR assistance

Real-time CPR assistance using metronome, compression velocity rendering and compression depth.

## Customisable

Voice messages, CPR settings, languages, defibrillation protocols, etc. can be customised via configuration tools.

## MobiMed Connect

Data transfer from the defibrillator to MobiMed Monitor and ePR.

## Self-tests

The defibrillator performs regular and automatic self-tests to ensure functionality and facilitate maintenance.

## Intuitive

Clear voice instructions and guiding icons via display.

## Patient history

The defibrillator has a memory function for storing patient cases.

## 360 - Joule

Energy levels up to 360 joules.

## Dust / water resistance

IP55 - Created for tough conditions



## Contact us!

Ortivus AB  
Svärdvägen 19  
182 33 Danderyd  
Sweden

[info@ortivus.com](mailto:info@ortivus.com)  
+46 0 446 45 00

