

## Release Notes MobiMed Smart v4.2.11

Software Article Number 490210837 (Software Build No 4.2.11\_151217.1)

### Customer specific changes

- SWASFT: CCN SW032 - Acute Authentication on CWS. (6029)
- SWASFT: CCN SW028: Saving Records as Other File Formats. Filename format changed to: <patient lastname>\_<patient forename>\_<NHS-Number>\_<CAD Incident Number>.<file extension>. (6143)
- SWASFT: Replaced non-time bound NEWS id's with time bound. (6205)

### Bugs corrected and improvements

- Bug corrected that could cause the program to hang when alerted patients where opened on a CWS. (6188)
- Closing expired TCP/IP sockets to reduce application server resource consumption. (6095)
- CAD information containing age sometimes caused an error message in the ePR, this has been corrected. (6194)
- Issue corrected which caused an error message to be displayed when the patient was discharged at the same time as the vital signs bar was updated. (6231)
- Issue with drawing pacemaker marks corrected, which could cause a crash. (6191)
- Issue addressed with having ' in group names for region, hospital and departments. (6160)
- Fixed inconsistency between DoB control and Age controls. (6195)
- It is now prepared for the setting of Stations on a PU to be connected to the settings tree, and hence have the stations listed instead of being entered as free texts in a free text field. The change simplifies for the user and prevents misspelled entries for Station. (5836)

### Known/Existing Limitations

- Handover between customers (i.e. between SCASFT and SWASFT) is not supported today and WiFi handover between different customers has been disabled. The side effect for this release is that it will not be possible to connect to WiFi from a PU with new version ( $\geq 4.2.11$ ) to an old version ( $< 4.2.10$ ). Vice versa works however. (6075)

*\*Numbers in brackets are ID-numbers from Ortivus internal source control, work item tracking, reporting and project management tool.*



## Required ePR version

The functionality in this release requires the below ePR version.

- SCASFT: SCASFT\_4.2.11\_2015.12.17.1
- SWASFT: SWASFT\_4.2.11\_2015.12.17.1

## Compatible hardware for ambulance units

The product is released for the following hardware platforms:

- Panasonic H2

All versions of Measuring Units M531 and 300C are supported.

## Release Information

The product is released on the following markets:

- UK – for below Ambulance Services
  - SWASFT (South Western Ambulance Services Foundation Trust)
  - SCASFT (South Central Ambulance Services Foundation Trust)

Replaces all previous 4.2 releases.

Upgrade to 4.2.11 can be done from all previous 4.2 releases.

This is a commercial release.