

## MobiMed Release Notes 4.53



Software part number: [490212514](#)

Software build number: 4.53\_2021.02.10.2

### Release Updates

The following are implemented in 4.53:

1. Change Requests implemented
2. General changes and improvements
3. Fixes



## Change Requests Implemented

### SWASFT SW2.43 - Dexamethasone [12560]

Updates to the medication list

### SCASFT SC312 - Hypoglycaemia [12134]

ePCR change to add templating for Hypoglycaemia incidents

### SCASFT SC317 - Stroke template [12135]

ePCR change to add templating for Stroke incidents

#### **SCASFT SC321 - Covid 19 and Oximeter [12529]**

ePCR change to enable registering of Oximeter devices left with a non-conveyed patient.

#### **SCASFT SC325 - Lateral Flow Test Tab [12530]**

ePCR change to register administered COVID-19 lateral flow test.

## **General changes and improvements**

#### **Remove obsolete version of TeamViewer [11282]**

Since support for Teamviewer has ended, a configuration has been added to remove the button to start Teamviewer from the user interface.

#### **Updates to the eprDateSelector [12149]**

Added possibility to set min and max times relative to the ePCR creation time.

#### **Integration engine updates [12181]**

Improvement of robustness and performance, by adding support for multi-threaded integration backends.

#### **Performance enhancement on DBI [12202]**


Improvement of performance by reusing network connections

#### **Ensuring that Ortivus Bluetooth stack works with Panasonic CF-20**

Adapting to handle Ortivus Bluetooth stack

#### **Improved review of ECG [11361]**

Added possibility to compare complexes.



## Fixes

Resolved text mismatch for label "Send 12" button in realtime view and in complex view [12152]

Resolved problem with swedish characters in data warehouse mission columns [12183]

Resolved issue with missing mandatory marks with Kronoberg's ePCR [12478]

Resolved issue with integration engine (HIS) regarding when deleting report on care episode that has multiple reports [12567]

Resolved issue with displaying etCO2 while using auto capture [12599]

Resolved issue in OSDK, error when configuration was missing [12600]

## Installation and deployment

For installing this release. Please contact our Service and Support at the below addresses.

- **UK market:** [support@ortivus.co.uk](mailto:support@ortivus.co.uk)
- **Other markets:** [support@ortivus.com](mailto:support@ortivus.com)

## Release Information

This is a commercial release on the following markets:

- Norway
- Sweden

- United Kingdom
- Italy

Upgrade to this release can be done from release 4.2 and onwards.

Please contact your customer representative for more information.

## Release Schedule



# MobiMed Release Plan

Release	Release date	Last feature request
4,52	21 Dec 2020	02 Nov 2020
4,53	15 Feb 2021	14 Dec 2020
4,54	29 Mar 2021	08 Feb 2021
4,55	10 May 2021	22 Mar 2021
4,56	21 Jun 2021	03 May 2021
4,57	30 Aug 2021	14 Jun 2021

### Ortivus

Svärdsvägen 19, Danderyd  
Sweden

+46 8 446 4500  
info@ortivus.com

You received this email because you signed up on our  
website or made a purchase from us.

[Unsubscribe](#)